



PIONEER RESEARCH CHEMICALS LTD
Manufacturer of Medical Reagents, Stains and Dyes



PRODUCT CATALOGUE

Issue: 08

Page

- 3 - H & E Stains
- 4 - Special Stains/Reagents
- 5 - Cytology Stains
- 6 - Stain Quick Solutions
- 7 - Pertex® Mounting Medium
- 8 - Fixatives & Wash Buffer
- 9 - Decalcifiers
- 10 - 10L Bag-in-a-box (various)
- 11 - Xylene Substitute
- 11 - Solvents

Page

- 12 - Davidson Marking Dyes
- 14 - Dry Stains & Dyes
- 15 - OCT Embedding Matrix
- 16 - Microbiology Stains
- 17 - Haematology Stains
- 18 - Reprotoxicology Stains
- 19 - Speci-wrap
(Histological Paper)
- 20 - Formalin Safety & Disposal
- 22 - Packaging




Pioneer Research Chemicals Ltd has established an excellent reputation, within the UK and internationally, for the manufacture and supply of high quality reagents, stains, fixatives and buffers. Our client base consists of hospitals, research facilities, distributors, universities, veterinary and general laboratories.


Established over 35 years PRC boasts a wide range of knowledge, expertise and experience. As an independent manufacturer we are able to offer complete flexibility, including the ability to supply custom made formulations. This flexibility enables us to respond quickly and effectively to customer requirements.


If you are interested in using or distributing our products please do not hesitate to contact us.

How to place an order or contact us:

 +44(0)1206 791781

 sales@pioneerresearch.co.uk

 www.pioneerresearch.co.uk

 7 Commerce Park, Whitehall Ind. Est., Colchester, Essex, CO2 8HX, UK


Newsletter Sign Up: Visit www.pioneerresearch.co.uk

Delivery Charges: Please see pricing table on the back page.

Terms and Conditions: Available at www.pioneerresearch.co.uk

Safety Data Sheets: Please email sales@pioneerresearch.co.uk. When requesting Safety Data Sheets please supply your name and company details.

HAEMATOXYLIN & EOSIN SOLUTIONS



EOSIN SOLUTIONS:

Product Code	Product	£ per Litre	£ per 2.5 Litre	£ per 5 Litre
PRC/66/1	Eosin (1% Aqueous)	11.80	25.00	46.00
PRC/R/68	Eosin (1% Alcoholic)	14.50	34.00*	62.00*

* **Note:** Hazardous delivery surcharge applies. See Pricing Table on back page for details.
Note: Other Eosin concentrations are available on request.

HAEMATOXYLIN SOLUTIONS:

Product Code	Product	£ per Litre	£ per 2.5 Litre	£ per 5 Litre
PRC/R/42	Haematoxylin (Mayers)	18.80	42.00	78.00
PRC/R/51	Haematoxylin (Harris) (Mercury Free) (Acidified)	19.20	44.00	82.00
PRC/R/51/2	Haematoxylin (Harris) (Mercury Free) (Unacidified)	19.20	44.00	82.00
PRC/13/2	Haematoxylin (Gills 1)	19.70	45.20	83.80
PRC/13/1	Haematoxylin (Gills 2)	19.70	45.20	83.80
PRC/13/3	Haematoxylin (Gills 3)	19.70	45.20	83.80
PRC/R/78	Haematoxylin (Coles)	18.40	42.30	79.00
PRC/R/131	Haematoxylin (Carazzis)	18.40	42.30	79.00
PRC/R/113	Haematoxylin (Weigerts) Iron Kit - contains 1L each of solution A and B	36.50		

SPECIAL STAINS & STAINING REAGENTS



SPECIAL STAINS:

Product Code	Product	Pack Size	Unit Price (£)
PRC/R/54	Van Gieson Stain	1litre	18.40
PRC/R/54/1	Van Gieson Stain (Curtis)	1litre	26.70
PRC/R/109	Picrosirius Red Stain	1litre	47.50
PRC/R/108	Millers Elastin Stain	1litre	58.70
PRC/16	Feulgen Stain (Schiiff's Reagent)	1litre	21.30
PRC/30	1% Alcian Blue in 3% Acetic Acid	1litre	48.00
PRC/6	1% Periodic Acid	1litre	16.00
PRC/R/135	0.2% Gold Chloride	500ml	120.00
PRC/25	Silver Nitrate 5% Solution	500ml	40.80
PRC/21	1% Neutral Red (Aqueous)	1litre	14.80
PRC/R/22	0.1% Nuclear Fast Red + 5% Aluminium Sulphate	100ml	16.60

FOR COLD STAINING OF ACID FAST BACILLI:

Product Code	Product	Pack Size	Unit Price (£)
PRC/41	Fuchsin Cold AFB Stain	1litre	23.80
PRC/42	Blue Counterstain Cold AFB Stain	1litre	17.70

STAINING REAGENTS:

Product Code	Product	£ per Litre	£ per 2.5 Litre	£ per 5 Litre
PRC/R/76	Scott's Tap Water Substitute	5.20	11.80	19.80
PRC/R/75	1% Acid Alcohol	8.20	18.80*	32.00*

* Note: Hazardous delivery surcharge applies. See Pricing Table on back page for details.



CYTOLOGY STAINS



PAPANICOLAOU STAINS:

Product Code	Product	£ per Litre	£ per 2.5 Litre	£ per 5 Litre
PRC/R/58	PRC Papanicolaou EA 50	15.40	35.00*	60.00*
PRC/R/77	PRC Papanicolaou OG 6	15.40	35.00*	60.00*

ROMANOWSKY STAINS:

Product Code	Product	£ per Litre	£ per 2.5 Litre	£ per 5 Litre
PRC/R/48	May Grunwald Stain	15.40	31.10*	55.00*
PRC/R/173	Giemsa Stain Modified	15.60	31.50	55.40

Giemsa Stain Modified: Contains no Methanol, safe to use, non toxic, stable, non flammable and no special storage required.

* **Note:** Hazardous delivery surcharge applies. See Pricing Table on back page for details.

STAIN QUICK (FOR RAPID ROMANOWSKY TYPE STAINING):


PRC Stain Quick is for rapid, reproducible routine staining of haematological and cytological material prepared onto microscope slides. **See page 6** for further Stain Quick information, including the suggested manual staining protocol.

PHOSPHATE BUFFER:

Product Code	Product	£ per 100ml
PRC/R/185	Phosphate Buffer (Sorensen) x50 Concentrate - makes 5 litres when diluted with de-ionised water. Available pH 6.8.	8.00

STAIN QUICK

FOR RAPID ROMANOWSKY TYPE STAINING



STAIN QUICK:

PRC Stain Quick is for rapid, reproducible routine staining of haematological and cytological material prepared onto microscope slides.

STAINING TECHNIQUE:

The technique can be completed in 2 minutes following the suggested manual staining protocol below:

- | | |
|---------------------------------|---------------------|
| 1 - Methanol | 10 dips x 1 seconds |
| 2 - Stain Quick Solution - red | 15 dips x 1 seconds |
| 3 - Stain Quick Solution - blue | 10 dips x 1 seconds |
| 4 - Buffer solution pH 6.8 | 2 dips x 10 seconds |

Further instructions are available at: www.pioneerresearch.co.uk/stain-quick/

STAIN QUICK KIT:

Product Code	Product	1 kit (£)	4 kits (£)
PRC/SQ/KIT	Stain Quick Kit: (100ml fixative, 100ml red solution, 100ml blue solution, 3 x pH 6.8 buffer tablets)	35.00	100.00

INDIVIDUAL SOLUTIONS:

Product Code	Product	Price per 500ml (£)	Price per 1 Litre (£)	Price per 5 Litre (£)
PRC/SQ/1	PRC Stain Quick Fixative	13.50	20.00	50.00*
PRC/SQ/2	PRC Stain Quick Red	13.50	20.00	50.00
PRC/SQ/3	PRC Stain Quick Blue	13.50	20.00	50.00

* Note: Hazardous delivery surcharge applies. See Pricing Table on back page for details.

PERTEX® MOUNTING MEDIUM



PERTEX® MOUNTING MEDIUM

Pertex is a mounting medium intended for use after clearing in xylene, ethanol and some xylene substitutes. Pertex is a quick-drying, xylene-based mounting medium with low viscosity. The low viscosity makes the medium flow out easily over the slide and prevents the formation of bubbles. A thin layer of Pertex is sufficient for mounting.

Pertex has the following attributes:

- Is invisible as the refractive index (1.492) is very close to the refractive index of glass
- Preserves the colour, form and structure of the slide over time
- Fast drying
- Compatible to xylene

Product Code	Product	Price per 250ml (£) (1 – 3 units)	Price per 250ml (£) (4+ units)
PRC/R/750	Pertex® Mounting Medium (250ml)	15.50	11.80

Stay connected with PRC – Sign up for our Newsletter at www.pioneerresearch.co.uk

Standing Orders can be set up so that pre-arranged quantities of products can be delivered at regular intervals.

FIXATIVES & WASH BUFFER

- Ready-made solutions saving time and money
- Custom formulations available
- Available in a range of pack sizes



FIXATIVES:

Product Code	Product	£ per Litre	£ per 5 Litre
PRC/R/1	Bouins Solution	9.80	32.00
PRC/R/4	10% Neutral Buffered Formal (pH 6.9 - 7.1)	9.80	14.90
PRC/R/43/C	Neutral Buffered Formal Concentrate 1:5	10.90	18.40
PRC/R/2	10% Buffered Formal Saline (pH 6.6 - 7.0)	9.90	15.50
PRC/R/69	10% Formal Saline (pH 3.5 - 4.5)	8.90	11.20
PRC/R/69/C	Formal Saline Concentrate 1:5	9.90	16.80
PRC/R/3	Davidson's Solution	9.90	32.00
PRC/R/200	Modified Davidson's Solution	12.00	38.50
PRC/R/850	Alcoholic Formalin	11.30	35.80*
PRC/R/505	Formaldehyde 38/10%	8.00	15.90
PRC/R/38	Parafix 4% (Electron Microscopy Fixative)	16.20	54.00
PRC/R/40	Karnovsky's Fixative (modified)	28.00	

WASH BUFFERS:

Product Code	Product	£ per Litre	£ per 5 Litre
PRC/R/89	Phosphate Buffered Saline (pH 7.2 - 7.4)	2.00	7.00

* Note: Hazardous delivery surcharge applies. See Pricing Table on back page for details.

All prices exclude carriage and VAT. See back cover for delivery charges.



DECALCIFERS



DECALCIFERS:

Product Code	Product	£ per Litre	£ per 5 Litre
PRC/R/12	D.F.B (Decalcifying Agent - H.K Kristensen)	6.80	23.00*
PRC/R/11	D.F.A (Decalcifying Agent - Gooding & Stewart)	9.10	38.00*
PRC/R/162	Formic Acid 10%	6.00	19.80*
PRC/R/162/1	Formic Acid 20%	6.60	21.40*
PRC/R/24	Formic Acid - Sodium Citrate Decalcifier	9.40	29.50*
PRC/R/118	EDTA	10.70	41.50
PRC/R/132	EDTA + 10% Formalin	12.00	45.00

* **Note:** Hazardous delivery surcharge applies. See Pricing Table on back page for details.

Bulk discounts on all Fixatives & Decalcifiers available

OEM Service – Choose PRC as your preferred outsourcing partner for private labelling & OEM services for histology solutions.

10L BAG IN A BOX

- **Reduced Environmental Impact**
- **PRC's 10L Bag-in-a-box uses 80% less plastic than our 10L jerry can!**



PRC's 10 litre bag-in-a-box pack provides users with:

- High Quality UN approved packaging for safe laboratory handling
- Ready-made solution for ease of use
- Dispensing Tap allows controlled dispensing
- Cube shaped for improved storage
- Ideal for bench top use
- Reduced environmental impact
- Reduced packaging disposal costs

FIXATIVES:

Product Code	Product Description	£ per 10 Litre
PRC/R/1/BB	Bouins Solution	50.00
PRC/R/69/BB	10% Formal Saline	19.60
PRC/R/2/BB	10% Buffered Formal Saline (pH 6.6 - 7.0)	24.40
PRC/R/4/BB	10% Neutral Buffered Formal (pH 6.9 - 7.1)	23.90

DECALCIFIERS:

Product Code	Product Description	£ per 10 Litre
PRC/R/12/BB	D.F.B (Decalcifying Agent - H.K Kristenson)	36.00*
PRC/R/118/BB	D.F.A (Decalcifying Agent - Gooding & Stewart)	57.00*
PRC/R/162/BB	Formic Acid 10%	31.80*
PRC/R/162/1/BB	Formic Acid 20%	34.50*
PRC/R/24/BB	Formic Acid - Sodium Citrate Decalcifier	42.40*
PRC/R/118/BB	EDTA	66.40
PRC/R/132/BB	EDTA + 10% Formalin	71.80

* Note: Hazardous delivery surcharge applies. See Pricing Table on back page for details.

WASH BUFFERS:

Product Code	Product Description	£ per 10 Litre
PRC/R/89/BB	Phosphate Buffered Saline (pH 7.2 - 7.4)	13.80

XYLENE SUBSTITUTE & SOLVENTS



XYLOPATH (XYLENE SUBSTITUTE):

Product Code	Product	Price per Litre (£)	Price per 5 Litre (£)
PRC/R/175	Xylopath	4.90	22.00

- Intended for use in histological laboratories as a xylene substitute
- Is safe to use, non aromatic, odourless, non toxic and non flammable. It can be used to reduce the health hazards associated with xylene
- Can be used with current automated tissue processing

Note: Xylopath is for histoprocessing only. It is not recommended for staining.

SOLVENTS:

Product Code	Product	Price per Litre (£)	Price per 5 Litre (£)
PRC/R/101	Industrial Denatured Alcohol (IDA) 99%	5.80	17.80*
PRC/R/102	Industrial Denatured Alcohol (IDA) 95%	5.80	17.80*
PRC/R/103	Industrial Denatured Alcohol (IDA) 70%	5.80	17.80*
PRC/R/201	Xylene (sulphur free)	5.40	18.90*
PRC/R/301	Toluene	8.00	21.50*
PRC/R/506	Isopropanol (IPA)	8.60	23.30*
PRC/R/680	Completely Denatured Alcohol (CDA)	5.90	18.50*

- CDA is a European approved formulation exempt from excise duties within the EU
- CDA users do not require HMRC User Authorisation
- CDA samples available for evaluation purposes

* **Note:** Hazardous delivery surcharge applies. See Pricing Table on back page for details.

Note: A copy of User Authorisation Letter is required prior to despatch of Industrial Denatured Alcohol products.

DAVIDSON MARKING DYES

Pioneer Research Chemicals Ltd. are pleased to have been appointed distributors of Davidson Marking Dyes for the UK and Ireland.



Davidson Marking Dyes are a permanent dye used to mark and orientate surgical and other tissue specimens. These high quality marking dyes enjoy a wide acceptance in hospitals, clinics and research facilities around the world.

Benefits include:

- Economical, efficient, visible, reliable and safe
- Permanently mark tissue for orienting and for other tissue marking applications
- Process multiple specimens in a single cassette
- Tissue Dye appropriate for both fresh and frozen tissue processing


Available in Green, Yellow, Black, Red, Blue, Orange and Violet

To view a video demonstrating the use of the Davidson Marking System search on You Tube for 'Using the Davidson Marking System'.

COMPLETE DAVIDSON MARKING SYSTEM:		
Product Code	Product Description	Price (£)
13-2406	Complete Davidson Marking System consisting of: 6 x 2oz bottles of dye (green, yellow, black, red, blue & orange) wooden holding tray and 50 applicator sticks.	170.00



DAVIDSON MARKING DYES



8oz bottles:

Product Code	Product	Price (£)
13-34081	8oz bottle green dye	74.00
13-34082	8oz bottle yellow dye	74.00
13-34083	8oz bottle black dye	74.00
13-34084	8oz bottle red dye	74.00
13-34085	8oz bottle blue dye	74.00
13-34086	8oz bottle orange dye	74.00
13-34087	8oz bottle violet dye	74.00

2oz bottles:

13-11631	2oz bottle green dye	37.00
13-11632	2oz bottle yellow dye	37.00
13-11633	2oz bottle black dye	37.00
13-11634	2oz bottle red dye	37.00
13-11635	2oz bottle blue dye	37.00
13-11636	2oz bottle orange dye	37.00
13-11637	2oz bottle violet dye	37.00

Miscellaneous:

13-44012	3cc 7 mini tube variety pack	10.00
13-8553	500 Applicator sticks (10 x packs of 50)	15.00
14-0126	Indian Ink (black) 500ml	46.80

DRY STAINS AND DYES



DRY STAINS AND DYES:

Product Description	CI Name	CI No. & Code No.	Price per 100g (£)	Price per 25g (£)
Alizarin Red S	Mordant Red 3	58005	119.80	34.30
Aniline Blue Water Sol.	Acid Blue 93	42780	53.20	20.70
Biebrich Scarlet	Acid Red 66	26905	30.20	9.10
Bismark Brown R	Basic Brown 4	21010	29.00	8.00
Eosin B	Acid Red 91	45400	35.90	10.30
Eosin Yellowish	Acid Red 87	45380	35.20	10.00
Fast Green FCF	Food Green 3	42053	66.00	22.00
Fuchsin Acid	Acid Violet 19	42685	47.50	14.40
Fuchsin Basic	Basic Violet 14	42510	30.80	10.30
Giemsa Stain	-	PRC/P/17	34.20	11.80
Haematoxylin	Natural Black 1	75290	138.00	38.50
Leishman Stain	-	PRC/P/18	34.70	11.70
Malachite Green	Basic Green 4	42000	34.70	11.70
May & Grunwald Stain	-	PRC/P/19	40.60	13.50
Methyl Violet 2B	Basic Violet 2B	42535	32.30	10.90
Methylene Blue BP73	Basic Blue 9	52015	30.00	9.30
Neutral Red	Basic Red 5	50040	53.00	15.60
Oil Red O	Solvent Red 27	26125	59.70	19.50
Orange G	Acid Orange 10	16230	29.70	9.80
Ponceau S	Acid Red 112	27195	48.00	16.00
Safranin O	Basic Red 2	50240	37.20	11.90
Sudan II	Solvent Orange 7	12140	31.50	9.00
Sudan IV	Solvent Red 24	26105	35.70	11.50
Victoria Blue B	Basic Blue 8	42563	27.70	8.90
Wrights Stain	-	PRC/P/28	72.00	24.30

Bulk discounts available - please contact sales@pioneerresearch.co.uk



OCT EMBEDDING MATRIX



Our OCT product provides a quality embedding matrix for frozen tissue sections.

Product advantages include:

- Rapid freezing
- Smooth cutting
- 125ml specially designed and easy to use stopper bottle
- Miscible with water, allowing easy removal
- Long shelf life
- Competitively priced

OCT EMBEDDING MATRIX:

Product Code	Product	Unit	£ per unit
PRC/OCT	OCT Embedding Matrix	125ml	10.50
PRC/OCT/4	OCT Embedding Matrix	4 x 125ml (1 box)	28.00

MICROBIOLOGY STAINS

(READY TO USE)



MICROBIOLOGY STAINS:			
Product Code	Product	Pack Size	Unit Price (£)
PRC/R/55	Carbol Fuchsin (ZN)	1 litre	15.90
PRC/41K	Carbol Fuchsin (Kinyoun)	500ml	17.70
PRC/41	Fuchsin Cold AFB Stain	1 litre	23.80
PRC/42	Blue Counterstain Cold AFB Stain	1 litre	17.70
PRC/R/137	Crystal Violet (0.5% Aqueous)	1 litre	14.80
PRC/R/53	Lugols Iodine	1 litre	15.30
PRC/R/203	Malachite Green (0.1% Aqueous)	1 litre	14.40
PRC/18	Methylene Blue (1% Aqueous)	1 litre	14.00
PRC/R/186	Methylene Blue Polychrome (MacFadyen)	1 litre	44.00
PRC/21	Neutral Red (1% Aqueous)	1 litre	14.80
PRC/44	Safranine (1% Aqueous)	1 litre	13.90

HAEMATOLOGY STAINS



HAEMATOLOGY STAINS

Product Code	Product	£ per Litre	£ per 2.5 Litre	£ per 5 Litre
PRC/R/48	May Grunwald Stain	15.40	31.10*	55.00*
PRC/R/173	Giemsa Stain Modified	15.60	31.50	55.40


Giemsa Stain Modified: Contains no Methanol, safe to use, non toxic, stable, non flammable & no special storage required.

* **Note:** Hazardous delivery surcharge applies. See Pricing Table on back page for details.

PHOSPHATE BUFFER:

Product Code	Product	£ per 100ml
PRC/R/185	Phosphate Buffer (Sorensen) x50 Concentrate - makes 5 litres when diluted with de-ionised water. Available pH 6.8.	8.00

REPROTOXICOLOGY REAGENTS



REPROTOXICOLOGY REAGENTS:

Product Code	Product	£ per Litre	£ per 5 Litre
PRC/R/25	Glycerol 99.5% minimum	9.60	38.80
PRC/R/25/1	Glycerol 75%	8.80	36.60
PRC/R/25/2	Glycerol 50%	8.00	32.80
PRC/R/25/4	Glycerol 25%	7.00	28.30
PRC/R/47	1% KOH in 20% Glycerol	6.20	28.00
PRC/R/67	1% KOH + 0.1% Alizarin (Kohal 8)	5.90	24.00
PRC/R/29	1% KOH (Potassium Hydroxide Solution)	4.60	20.00
PRC/R/19	2% KOH (Potassium Hydroxide Solution)	5.10	22.00

Other formulations are available to suit your requirements, for further details please contact us.

Custom Solutions – our in-house knowledge and flexibility means we are able to deliver custom made formulations.

SPECI-WRAP (HISTOLOGICAL PAPER)

Speci-wrap for wrapping small samples and biopsies during processing.



SPECI-WRAP:

Product Code	Product	£ per box
11-WR60	Speci-wrap - 60 x 60mm (Box of 500 sheets)	9.60

Stay connected with PRC – Sign up for our Newsletter at www.pioneerresearch.co.uk

Pricing – whilst we endeavour to maintain our prices throughout the life of this brochure, we reserve the right to change them without notice. All prices exclude VAT and delivery charges. If you are VAT exempt, please include a copy of your exemption certificate with your purchase order.

FORMALIN SAFETY & DISPOSAL



FORMALDEHYDE NEUTRALISER PADS:

Designed for lining trays, cassette inspection and examination and dissection of biological specimens. The pad quickly absorbs and neutralises formaldehyde. Available in 4 convenient pack sizes.

Product Code	Product	Price (£)
PRC/FNP/S	Formaldehyde Neutraliser Pad - 20 x 25 cm (8" x 10")	33.90 (25 pads)
PRC/FNP/M	Formaldehyde Neutraliser Pad - 28 x 35cm (11" x 14")	59.50 (25 pads)
PRC/FNP/L	Formaldehyde Neturaliser Pad - 40 x 50cm (16" x 20")	102.80 (25 pads)
PRC/FNP/R	Formaldehyde Neutraliser Roll - L: 18m x W: 34cm	134.00 (each)

FORMALIZE (WASTE FORMALIN NEUTRALIZER):

Formalize neutralizes 4% formaldehyde solutions overnight, by converting them into non-hazardous products. Supplied in buckets, sachets and ready-to-fill 10 litre containers, Formalize improves safety and helps with waste disposal management.

Product Code	Product	Price Each (£)	Price per Pack (£)
PRC/R/14/B	Formalize bucket (to neutralize 10 litres of 4% formaldehyde)	6.60	29.70 for 6 (4.95 each)
PRC/R/14/S	Formalize sachet (to neutralize 5 litres of 4% formaldehyde)	4.40	82.00 for 25 (3.28 each)
PRC/R/14/C	Formalize (ready-to-fill 10 litre container + tap)	9.80	N/A



FORMALIN SAFETY



SPILL KIT (4% FORMALDEHYDE SPILL):

Response kit for absorbing and neutralising 4% formaldehyde spills which may occur in transit and laboratory areas.

The kit contains:

- 1 tub Formaspill Absorbent (500g)
- 1 mask - disposable
- 1 pair safety glasses – disposable
- 2 pairs nitrile gloves – disposable
- 1 plastic apron – disposable
- 1 pair plastic sleeves – disposable
- 1 plastic scoop and scraper – disposable

Product Code	Product	Price (£)
PRC/FSK	4% Formaldehyde Spill Kit	28.00

FORMASPILL ABSORBENT:

Formaspill rapidly absorbs and neutralizes 4% formaldehyde, allowing spills to be safely contained, with the elimination of fumes and vapours.

Product Code	Product	Price Each (£)	Price per Pack (£)
PRC/R/95	Formaspill Absorbent (500g)	10.40	99.00 for 12 (8.25 each)

PACKAGING



CONTAINERS:

Product Code	Container Description	£ per Container
PRC/CONT/150	150ml cylindrical bottle, plastic, natural, 28mm Neck plus plastic, natural cap spouted with strap. Ideal for dispensing liquids.	7 for 10.00
PRC/CONT/1UN	1 Litre, UN, fluorinated, round, natural plus 42mm T/E wadded cap	1.20
PRC/CONT/5	5 Litre Jerry Can, natural plus 38mm wadded cap	1.05
PRC/CONT/5UN	5 Litre Jerry Can, UN, fluorinated, natural, plus 42mm T/E wadded cap	2.40
PRC/CONT/10UN	10 Litre Jerry Can, UN, natural, stackable plus 51mm T/E cap	3.55
PRC/CONT/25UN	25 Litre Jerry Can, UN, blue, stackable plus 61mm T/E cap	7.10

WIDE NECK, RECTANGULAR, DOUBLE CLOSURE SCREW CAP BOTTLES

- Heavy duty, rectangular, wide neck, HDPE natural bottle with plug, seal and screw cap
- Easy to grip caps
- Moulded graduations for approximate volume

Product Code	Capacity (ml)	Neck (mm)	L x W x H (mm)	£ per bottle
PRC/WN/500	500	45	80 x 65 x 135	2.10
PRC/WN/1000	1000	58	100 x 80 x 177	3.20
PRC/WN/2000	2000	58	132 x 100 x 214	5.75

Labelling facility available on request. Price on application.

CARTONS:

Product Code	Description	Internal (LxWxH) mm	Pack Qty	£ per pack
PRC/B/1	Single wall to fit 4 x 1 Litre containers	177 x 177 x 241	10	10.00
PRC/B/3	Double wall to fit 4 x 5 Litre containers	380 x 285 x 295	15	21.00

All prices exclude carriage and VAT. See back cover for delivery charges.



DELIVERY CHARGES

	DESTINATION		
	England *, Wales & Scotland (exc. VAT)	N.Ireland, Rep. of Ireland, Channel Islands & Isle of Man (exc. VAT)	Other EU countries
Non Hazardous - Up to 20 Litres	£10.00	£37.50	POA
Hazardous - Up to 20 Litres	£35.00	POA	POA
Non Hazardous - 21 to 40 Litres	£17.50	POA	POA
Hazardous - 21 to 40 Litres	£55.00	POA	POA
Non Hazardous - 41 to 60 Litres	£25.00	POA	POA
Hazardous - 41 to 60 Litres	£65.00	POA	POA
Non Hazardous - 61 to 500 Litres	£45.00	POA	POA
Hazardous - 61 to 500 Litres	£70.00	POA	POA

- All consignments 500+ litres are price on application (POA)
- * Deliveries to the Isle of Wight will incur a £6 delivery surcharge.
- We only use reputable national courier companies that have track and trace services.
- If you have any questions relating to our delivery charges please contact us.

Prices

All prices exclude VAT and delivery charges. If you are VAT exempt, please include a copy of your exemption certificate with your purchase order.

Whilst we endeavour to maintain our prices throughout the life of this brochure, we reserve the right to change them without notice.

